

How to Setup iRVision

1. Login as an Administrator
2. Copy the folder **iRVision Controls for R-30 iA V8.xx** from *T://IT/SOFTWARE-ISO/Niles - Robotics lab - iRVision* to **C:\Users\Public\Public Downloads** folder.
3. Logout.
4. Login as SMCSTUDENT.
5. Do not have any other programs running.
6. Expand the **vsctrls** zip file located in **the C:\Users\Public\Public Downloads\Niles - Robotics lab - iRVision**.
7. Run the **vsctrls** installation file in the **C:\Users\Public\Public Downloads\Niles - Robotics lab - iRVision\vsctrls** folder as Administrator.
8. When prompted, choose **Reinstall Using Recommended Settings**.
9. When the Security Warning shows, click on **RUN**.
10. Start the **HandlingPRO** program by clicking on it on the Desktop.
11. When it has started, click on **OPEN CELL** or click on the open folder.
12. In the folder **C:\Users\smcstudent\My Documents\My Workcells\HandlingPRO1-calibration** click on the file **HandlingPRO 1**
13. Unplug the Internet, and plug in the camera.
14. Start Internet Explorer (NOTE: Other browsers **will not work**).
15. Click on the gear at the top right-hand corner and choose **Compatibility View Settings**.
16. Add **localhost** to the list if it is not there already.
17. **CLOSE** the Compatibility Windowo.
18. Type **localhost** in the URL bar and hit **ENTER**.
19. Click on **Vision Setup** under the **iRVision** section.
20. If the Windows Firewall shows up, click on **ALLOW ACCESS**.
21. Double-click on **Test Camera**.
22. Ensure that the **Serial Number** on the right-hand side shows a number. If not, click the down arrow and select it.
23. You should see a picture from the camera. To see a **live** photo, click the small green screen icon at the top left under the menu. (Click on the red one to stop).
24. Click on **SAVE** at the bottom of the screen.